

Rultract Inc., 5663 Brecksville Road, Cleveland, Ohio 44131, U.S.A Website: [www.rultract.net](http://www.rultract.net)  
 East Coast Phone (440) 255-2910, Fax: (440) 255-2923 West Coast Phone (949) 725-9320, Fax: (949) 725-9751



**Rultract® Interchangeable Rakes**

**Pediatric Rakes**



**4100-60AS**  
 Pediatric Rake  
 Short, with 3 sharp tips  
 Full radius, 3/32 x 3/4  
 Stem 2 1/4" long  
 (2 mm x 19 mm)

**4100-60P-01**



Pediatric Rake  
 Short Stem, with 3 blunt tips  
 Full radius, 3/32 x 3/4  
 Stem 2 1/4" long  
 (2 mm x 19 mm)

**4100-60P-04**



Pediatric Rake  
 Single prong, blunt tip  
 Full radius, 1/8x 1/2  
 Stem 4" long  
 (3 mm x 13 mm)

**4100-60P-02**



Pediatric Rake  
 Short Stem, with 3 sharp tips  
 No radius, 3/32x 3/4  
 Stem 2 1/4" long  
 (2 mm x 19 mm)

**4100-60P-05**



Pediatric Rake  
 Single prong, blunt tip  
 Right angle bend, 1/8x 3/8  
 Stem 4" long  
 (3 mm x 9.5 mm)

**4100-60P-03**



Pediatric Rake  
 Short Stem, with 2 blunt tips  
 Full radius, 1/8x 1/2  
 Stem 2 1/4" long  
 (3 mm x 13 mm)

**4100-60P-06**



Pediatric Rake  
 Single prong, blunt tip  
 Full Radius, 1/8x 5/16  
 Stem 4 5/8" long  
 (3 mm x 8 mm)

